

Pyle Health
Bluetooth 4.0
Infrared Ear Thermometer



MODEL: PHTM10BT
User Manual

Contents

Cautions	2
Special Features	2
Measurement Principle	2
Product Application	3
Product Description	3
Ear Measurement	4
Connect device via Bluetooth 4.0	5
Correct Measurement Method	6
When the Measurement Result is Abnormal	6
Memory Function	7
Replacing the Batteries	7
Cleaning and Care	8
Product Specifications	8
Display Message	9

Thank you for purchasing infrared ear thermometer. Please read through this manual before using the product.

Cautions

1. Do not drop or damage the unit.
2. Do not immerse the unit into water.
3. Keep the unit out of reach of children.
4. If the unit does not function, contact the manufacturer.
5. Consult doctor if the result of measurement is out of normal body temperature range.
6. Do not disassemble the unit. Open the battery cover only when replacing battery.
7. Dispose used battery according to local recycling regulation.
8. The product function maybe reset when it is subjected to electrostatic discharge disturbance test. Just follow the instruction described in this manual to set the function again for normal operation.


Special Features

1. One second measure time.
2. LCD display for easy measurement reading.
3. Backlight display.
4. Memory for 30 readings.
5. Auto power off.

Measurement Principle

Because the eardrum shares blood supply with hypothalamus which is the temperature control center in the brain, measuring the temperature of the eardrum with infrared

thermometer can accurately reflect the body temperature.




Normal body temperature range according to World Health Organization(WHO)

ear	35.8°C - 38°C
Axillary	34.7°C - 37.3°C
Oral	35.5°C - 37.5°C
Rectal	36.6°C - 38°C

Product Application

For measurement of body temperature through eardrum.

Product Description



1. Probe
2. LCD display
3. Measurement button
4. Power/memory button
5. Battery cover

Ear Measurement

1. Press Power/Memory button Ready signal beep + All display + Backlight + last measurement result + Ready display The "°C" symbol is blinking.
2. Start measuring Pull ear smoothly to flatten ear canal and place the probe facing eardrum in straight line. Press measurement button until you hear the beep signal.


How to pull the ear




- For babies below 12 months: Pull ear smoothly backward.
 - For children above 12 months old and adults: Pull ear smoothly upward.
3. Measurement result When the measurement is complete, the temperature will be displayed with backlight.
 4. Interval Between Measurements If you need to take another measurement, allow an interval of 8 seconds until the signal "°C (°F)" flashes. If you press the Measurement button within 8 seconds, the beep signal will be sounded.
 5. Turn off the thermometer
 - Push Power/memory button until the "OFF" signal is displayed.
 - The thermometer will turn off automatically without use for 1 minute.

Connect device via Bluetooth 4.0

1. Installation Prior to first use, Download and install "Pyle Health" application to your device. (Bluetooth 4.0 capabilities, e.g. iPhone 4S, iPhone 5, NewiPad, iPadMini or Galaxy S4 etc..)
2. Open "Pyle Health" application.
3. Press the "Enter" button.




4. Press the "Thermometer" icon.



5. Press the "U" to turn on Thermometer When the measurement is complete, the temperature will display on App and save on iOS application screen.




6. The result will stores the last temperature reading:




Correct Measurement Method

1. Ensure the probe is facing eardrum in straight line. Otherwise, the measured temperature will be lower than the actual body temperature.

36.8°C



Lower measured temperature



2. For accurate result, clean the probe with cotton bud before measuring.


When the Measurement Result is Abnormal

1. Possible causes for high temperature:
 - After exercise, meal, or crying.
 - Prolonged exposure to high temperature environment.
 - Thermometer is taken out from low temperature environment. Keep the thermometer in room temperature for at least 30 minutes before using.
2. Possible causes for low temperature:
 - The probe is dirty.
 - The thermometer is pulled out from the ear before measurement complete.
 - The probe is not fitted snugly into the ear canal.
 - Just exposed to cold temperature environment before taking measurement.

- The probe is not facing the ear drum or is not placed in a correct position in the canal. See "Correct Measurement Method".
- Excess ear wax build-up.
- 3. Temperatures measured from left and right ears are different: There could be slight difference in temperatures when taking from both ears. Therefore, always take the temperature from the same ear.

Memory function

This thermometer stores the last 30 temperature readings:




- When powering on, press power/memory button to recall first reading.
- Repeat pressing power/memory button to recall the next reading.
- Delete memory: Press both power/memory and measurement buttons and hold for 3 seconds. The signal "Clr" displayed and the memory will be deleted.



Replacing the Batteries

This thermometer is powered by 2 AAA batteries. When the power is low, the low battery signal is displayed. To replace batteries: 1. Press and push the battery cover to open it.

2. Remove the old batteries and replace with new ones. Make sure the batteries are in right polarity direction.
3. Close the battery cover.



Note:

To protect the environment, dispose of empty batteries at your retail store or at appropriate collection sites according to national or local regulations.

Cleaning and care

1. Use a soft, dry cloth to clean the exterior of the thermometer.
2. Keep the thermometer away from direct sunlight, dust and contamination environment.
3. Store the thermometer in normal room temperature, or within 10-40°C.

Product Specifications

Measurement method	Infrared measurement
Measurement range	34-42.2°C
Accuracy	+/-0.2°C (ambient temp. 16-35°C) +/-0.3°C (ambient temp. 10-16°C, 35-40°C)
Measuring time	1 second
Display	LCD display
Memory	Stores last 30 readings
Auto power off	1 minute after not being used
Batteries	2 x AAA batteries
Operating temperature	10-40°C
Storage temperature	-20~+55°C
Dimensions	5.8" x 1.4" x 1.2" inch
Weight	1.73 oz. (not including batteries)

Display Message

Display message	Situation	Solution
	Ambient temperature is out of operating range (10-40 °C)	Allow the thermometer remain in room temperature (10-40 °C) for 30 minutes.
	Sensor error	Contact the dealer.
	Memory error	Contact the dealer.
	Temperature measured is higher than normal human temperature range (above 42.2 °C)	Check the measured object and take a new measure.
	Temperature measured is lower than normal human temperature range (below 34 °C)	Check the measured object and take a new measure.
	Low battery	Replace batteries.