

# Hot Stone Massage Kit

## Portable Heated Rock Massaging Therapy System

SereneLife

**USER GUIDE**

**PLEASE MAKE SURE THAT YOU HAVE READ THE MANUAL CAREFULLY BEFORE USING IT**

### **Hot Stone Massage Kit**

Thermal therapy, is the use of heat in therapy for pain relief and health. It can take the form of a hot stone, hot bamboo sticks, and others. It can be beneficial to those with arthritis and stiff muscles and injuries to the deep tissue of the skin. Our innovative electric massage warmer set comes with a lightweight energy efficient stylish storage/ travel bag case that doubles as stone heater oven that quickly heats up the rocks. It is a helpful tool to heat the massage stones and also bamboo sticks, as we know the application of warmed Basalt stones/bamboo sticks combined with manual massage to provide the client with a deeply relaxing, therapeutic experience. This treatment which also benefits the therapist not only financially but also physically.

### **Advantages of Hot Stone Massage Kit**

**Easy to Use:** The entire system is very easy to install and does not require water to heat the stone. Simply put the stone into the box and just plug in and set your temperature and time.

**Innovative Waterless Heating Technology:** With waterless heating technology there is no need for traditional stone heating equipment: thermometers, gloves, spoons, nets, etc. Due to the accumulation of scale, there is no need for regular cleaning. The stone is ready to be used about 15 minutes. Unlike stone heaters that use water, there is no mess after use.

**Portable Massage Kit:** Our portable massage kit offers mess-free stone warming along with the ability to easily transport for traveling therapist as well as your massage practice or home use. Our lighter weight heating case can heat stones and be transferred while retaining heat to make things easier on the therapist/user.

**Increased Productivity:** Not only does the Bag take less time to set up, but unlike a water stone heater, the bag can be transported whilst retaining heat even after its been disconnected from the mains, allowing the therapist to preheat the stones at home, thus saving time (and money!) at the clients house. See our quick formula (right) for how the bag pays for itself.

**Digital Temperature Controller:** Optional temperature setting is 85-160°F.

You just plug it in and you have hot stones, ready to be used at the right temperature. This feature is unique to Hot Stone Massage Kit.

**Multi functionality:** Combined heating with bag, it is more than a bag, not only can self-contained but also can heat.

**Heated Suitcase:** The electric massage heater set is equipped with a storage/travel bag that doubles as a stone heating furnace to quickly heat rocks. The massage warmer is the perfect choice for home use and travel, and is the perfect complement to travel massage therapists.

### How to use:

#### FOR 110V

**STEP 1:** Place the stones inside the heating warmer.

**STEP 2:** Connect the temperature controller to the the heating warmer

**STEP 3:** Connect the power cord to the temperature controller.

**STEP 4:** Press "ON/OFF" then can set the temperature from 85-160°F

#### "ERO"

1. Temperature controller has been connected power but not the heating warmer, in this situation, connect the heating warmer then it will work properly.
2. Temperature controller has been connected power as well as the heating warmer and NTC disconnected.

**GREEN LIGHT ON:** Heating system is working

**RED LIGHT ON:** Power has been connected

#### WARNING

1. Should be stored in a clean and dry environment. Wipe the Hot Stone Massage Kit with a damp cloth. To protect against electrical hazard, do not immerse in water or other liquid.
2. This product requires professional use.  
Please do not separate, maintain, or refit the product if you were not a professional.
3. The Kit should be keep away from high temperature environment and point-blank sunshine.
4. Do not leave unattended
5. Keep away from children!
6. Risk of electric shock, can cause injury or death if improper operation.
7. Do not block any ventilation when using, it will cause the fire in extreme cases.

#### Questions? Issues?

We are here to help!

Phone: (1) 718-535-1800

Email: [support@pyleusa.com](mailto:support@pyleusa.com)

#### California Prop 65 Warning

##### WARNING:

This product contains Wood Dust which is known to the state of California to cause cancer birth defects and other reproductive harm. Do not ingest.

For more info go to: [www.P65warnings.ca.gov](http://www.P65warnings.ca.gov)