

SereneLife
by **PYLE**

SereneLife
by **PYLE**



VISIT US ONLINE:

Have a question?

Need service or repair?

Want to leave a comment?

PyleUSA.com/ContactUs

PYLE



SLFTMSG18

Shiatsu Foot Massager

Deep Kneading Foot Massage with Heat

CONTENTS

- PRODUCT PARTS
- USING METHOD
- SAFE ATTENTION
- FEATURES
- TECHNICAL SPECIFICATIONS

Please read the instructions carefully before operating this machine, and keep it for later use.

NOTE: We reserve the right to revise or explain the design of this product without further notice.


• PRODUCT PARTS

A. Power Switch


C. Heat Lamp

B. Massage Heads

D. Cloth Cover



• REMOTE CONTROL FUNCTION



• USING METHOD

- Plug power cord to a 110-220V Ac outlet - massager will be in standby mode.
- Press the power switch, massage modes as follows:
 - **1st Press:** Strong kneading with heat
 - **2nd Press:** Slow kneading with heat
 - **3rd Press:** Strong kneading without heat
 - **4th Press:** Slow kneading without heat
 - **5th Press:** All functions off (unit returns to standby mode)

NOTE:

- A 15 minute massage is recommended. The unit will automatically stop after 15 minutes.
- The massage may tickle, which is normal that means sign of blood circulation improvement.
- Adjust the placement of your feet on the foot beds to massage different points of the feet.
- **IMPORTANT:** Do not sit or stand on the unit or exert excessive pressure on the foot beds. This may cause overheating and damage to the motor.

- **SAFE ATTENTION**

Usage Environment

- Do not use the massager in an excessively humid environment such as the bathroom.
- Do not use the massager on a wet body, This may result in short circuit or electric shock.

Unsuitability of User

The following user must consult their doctor before use

- Those suffering from cancer.
- Those suffering from osteoporosis.
- Those suffering from malignant tumors.
- Those who have a high fever.
- Women who are pregnant or menstruating.
- Persons who are insensitive to heat.

Safety Precautions during Use

- Do not use the massager for any purpose other than massage.
This may result in accident or injury.
- Do not use massager when foot pad fabric cover is torn or parts are missing

- Do not let children use without adult supervision or persons not capable of handling and operating the massager. Keep the massager out of reach of children.
- Do not stand or sit on the massager.
- Do not insert metallic objects, or flammable materials into the main body. This may result in electrocution, shot circuit and possible fire.
- Do not use the massager within an hour of meal time.
- Do not cover the massager during operation. This may result in a rise in motor temperature and subsequent breakdown or fire.
- Do not unplug or switch off the massager during operation. This may result in injury.
- Do not sleep while using the massager.

Assembly and Repair

- The massager must only be serviced or repaired by the Service Center. Users must not disassemble or repair the massager by themselves.

Instances when the massager shouldn't be used

In the following events, please switch off the power immediately and disconnect the AC power cord from the electrical outlet.

- If water is accidentally spilled on the massager.
This may result in electrical shock.
- If severe or abnormal pain is felt during massage, stop operation immediately and consult your doctor.
- If you detect any malfunction or any other abnormal condition during operation.
- If there is a power failure, injury may occur when power is restored unexpectedly.
- If there is lightning.

Notes about Power Plug and Cord

- Check that the voltage corresponds with specification indicated for the massager.
- Do not disconnect the AC plug from the electrical outlet roughly.
- Do not connect or disconnect the AC plug from the electrical outlet with wet hands. It may result in electrical shock or machine trouble.
- When disconnecting the AC plug, pull it out by holding the plug and not the cord.
- Do not place the AC cord under the massager or any heavy object.
- Do not wind the AC cord around the massager as it may damage the cord and result in fire or electrical shock.
- Do not operate the massager with a damaged cord or tamper with the cord. Contact the Service Centre for inspection.
- Please use the adapter and attachment we recommend.
- Do not use the massager if the electrical outlet is excessively loose.

Product Care and Storage

- Be sure to remove the AC plug from the electrical outlet before cleaning the massager.
- Do not use benzene, thinner or any other abrasive solvent on your massager.
- Use a dry cloth when cleaning the massager. Never use damp cloth to clean the massager or the power switch.
- When necessary, detach the fabric covers for cleaning to maintain good hygiene
- When not in use for a long period of time, it is advisable that the AC power cord of the massager be coiled and the massager should be stored in a dust and moisture free environment.
- Do not store the massager near intense heat or an open flame
- Do not leave massager under direct sunlight for long period of time.

SLFTMSG18

Shiatsu Foot Massager Deep Kneading Foot Massage with Heat

Features:

- Soothing Massage & Heat, for Your Feet!
- Selectable Deep Kneading Settings
- Increase Blood Flow & Circulation
- Relaxes & Relieves Muscle Pain
- Simple Electric Plug-in Operation
- Button Activated Heat Therapy
- Used for Toes, Feet, Heel & Calf Massaging
- High-Powered Oscillating Penetrating Motors
- Control Massage Settings via Remote Control

What's in the Box:

- Foot Massager
- Remote Control
- Power Adapter

Technical Specs:

- Heat Ability: ON / OFF
- Auto Timer: 15 mins / ON / OFF
- Selectable Kneading Settings: High / Low
- Power: 110/220V (24V Power Adapter)